

Applicable to

- Looking for a natural lifting
- Those who want skin care without interrupting daily life
- Concerned about decreased skin elasticity
- Wanting more meticulous skin care

Cartridge Types

	Cartridge	Area
FACE	1.5mm	Eye rim, forehead
	2.0mm	Cheekbone, cheek
	3.0mm	Double chin, neck, cheek
	4.5mm	Cheek
BODY	7.0mm	Arm
	10.0mm	Thigh
	13.0mm	Abdominal



ONE THERA

This product is not a medical device.
It is a reference material to help understand the product.

ONE THERA

DUAL MODE ULTRASOUND SYSTEM


What is One Thera?

One Thera is a device that delivers concentrated thermal energy to the dermis and fat layer of the skin through various types of cartridges. It can be used for both FACE & BODY, and it is recommended to select a cartridge and mode suitable for the area to be used and operate in an optimal way.



Long Pulse Mode?

Technology that divides one pulse into multiple pulses to deliver deep heat and further reduce pain.




● ONE THERA ● OTHERS

Normal mode & Long pulse mode

Normal: Delivers thermal energy to the desired skin layer
 Long pulse mode: Transmits thermal energy deeper

- Approximately 4 minutes for 3,000 shots up to 12Hz speed
- Various applications possible by using two modes



01 When can I see the results?

Normal: Delivers thermal energy to the desired skin layer
 Long pulse mode: Transmits thermal energy deeper

02 How often can I receive it?

Depending on the individual, you can expect immediate results immediately after use, and you will see better results after 2 to 4 weeks.

03 How is the pain level?

Depends on the individual, but generally it is a tingling degree of pain. If you are sensitive to pain, lower the energy level.

04 What are the precautions?

If the skin is red, it usually subsides on the day of treatment. It is recommended to refrain from sauna, drinking, and smoking for 3 to 5 days after the procedure.